


# WOUND HEALING PHASES


## A Practical Guide



### T TISSUE MANAGEMENT (CLEANING)

Procedure:

- Cleansing of wound bed
- Hydration of the wound
- Mechanical or surgical debridement
- Moisture management



Wound cleansing and wound bed preparation

Surgical / mech. Debridement, wound irrigation solution, e.B. ProntoVet® Prontosan® Debridement Pad

Primary dressing

Wound gel e.B. ProntoVet® Gel oder Prontosan® Gel X

Secondary wound dressing


Bandage e.g. Askina® Haft Color

### I INFLAMMATION, INFECTION

Procedure:

- Wound cleansing
- Antimicrobial cleaning
- Odor control
- Frequency of change of dressing: daily

Note: If there is no infection normal wound pads such as Askina® Foam/ Sorb can also be used.



Wound cleansing and wound bed preparation

Wound irrigation solution, e.g. ProntoVet®

Primary dressing

(Antimicrobial) wound dressing e.g. Askina® Calgitrol®, Askina® Foam, Askina® Sorb

Secondary wound dressing

Bandage e.g. Askina® Haft Color

## G GRANULATION (HUMIDITY CONTROL AND BALANCE)

Procedure:

- Mechanical / surgical debridement when necessary
- Wound cleansing
- Exudate management
- Granulation enhancement
- Frequency of change of dressing: approx. 1-3 days
- Maintain moisture balance  
(moisten dry wound additionally using ProntoVet® Gel, for instance)



EARLY GRANULATION



LATE GRANULATION

Wound cleansing and wound bed preparation

Surgical / mech. Debridement, wound irrigation solution, e.B. ProntoVet®

Primary dressing

Absorbing dressing (alginate / foam, foam / silicone) e.g. Askina® Sorb, Askina® Foam, later: Askina® DresSil

Secondary wound dressing

Bandage e.g. Askina® Haft Color

## E R EPITHELIALIZATION, REPAIR

Procedure:

- Wound cleansing
- Maintain moisture balance
- Protection of newly formed tissue
- Frequency of change of dressing: approx. 1-5 days



Wound cleansing and wound bed preparation

Wound irrigation solution, e.g. ProntoVet®

Primary dressing

Low adherence (silicone) dressing, hydrocolloid dressing (polyurethane film) e.g. Askina® DresSil / SilNet Plus

Secondary wound dressing

Bandage e.g. Askina® Haft Color