
Easyump® II ST/LT

en	Instruction for use	10
de	Gebrauchsanweisung	12
bg	Инструкция за употреба	14
cs	Návod k použití	16
da	Brugsanvisning	18
el	Οδηγίες χρήσης	20
es	Instrucciones de uso	22
et	Kasutusjuhised	24
fi	Käyttöohje	26
fr	Mode d'emploi	28
hu	Használati utasítás	30
id	Petunjuk penggunaan	32
it	Istruzioni per l'uso	34
kk	Қолдану нұсқаулығы	36
lt	Naudojimo instrukcija	38
lv	Lietošanas instrukcija	40
ms	Arahan Penggunaan	42
nl	Gebruikersinformatie	44
no	Bruksanvisning	46
pl	Instrukcja użytkowania	48
pt	Instruções de utilização	50
ro	Instrucțiuni de utilizare	52
ru	Инструкция по применению	54
sk	Návod na použitie	56
sl	Navodila za uporabo	58
sv	Brugsanvisning	60
tr	Kullanım Kılavuzu	62
uk	Інструкції для застосування	64
vn	Hướng dẫn sử dụng	66
zh	使用说明	68

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OVER- AND UNDERFILLING OF EASYPUMP® II ST-LT

PRODUCT SPECIFICATIONS

ST-VARIANTS										
ST-Variants	100 - 0.5	250 - 0.5	50 - 1	100 - 1	250 - 1	400 - 2	500 - 2	100 - 2	400 - 4	
Code No. (REF)	4540040	4540042	4540044	4540046	4540048	4540050	4540052	4540054	4540056	4540058
Nominal Flow Rate (ml/h)	200	500	50	100	250	175	200	250	50	100
Nominal Fill Volume (ml)	100	250	50	100	250	250	400	500	100	400
Minimal Fill Volume (ml)	50	125	30	125	135	135	240	240	50	240
Maximal Fill Volume (ml)	125	295	65	125	335	335	560	560	125	560
Maximum Retained Volume (ml)	≤3	≤8	≤2	≤3	≤8	≤8	≤10	≤10	≤3	≤10

STORAGE

TIME TO REACH ROOM TEMPERATURE (APPROX. +23 °C ±2 °C)										
Hours From Refrigerator (approx. +2 °C to +8 °C)	6	12	6	6	12	12	12	12	6	12
Hours From Freezer (approx. -18 °C)	12	18	12	12	18	18	18	18	12	18

OVER- AND UNDERFILLING DATA

APPROXIMATE DELIVERY TIME										
Hours (hh:mm)	Fill Volume (ml)									
00:15 h	50	125	-	-	-	-	-	-	-	-
00:30 h	100	250	-	50	-	-	-	-	-	-
00:45 h	-	-	38	75	188	-	-	-	-	-
01:00 h	-	-	50	100	250	175	-	250	50	-
01:15 h	-	-	63	125	-	219	-	313	63	-
01:30 h	-	-	-	-	-	263	300	375	75	-
01:45 h	-	-	-	-	-	-	350	438	88	-
02:00 h	-	-	-	-	-	-	400	500	100	-
02:15 h	-	-	-	-	-	-	450	-	113	-
02:30 h	-	-	-	-	-	-	500	-	125	250
02:45 h	-	-	-	-	-	-	550	-	275	-
03:00 h	-	-	-	-	-	-	-	-	-	300
03:30 h	-	-	-	-	-	-	-	-	-	350
04:00 h	-	-	-	-	-	-	-	-	-	400
04:30 h	-	-	-	-	-	-	-	-	-	450
05:00 h	-	-	-	-	-	-	-	-	-	500
05:30 h	-	-	-	-	-	-	-	-	-	550

PRODUCT SPECIFICATIONS

LT-VARIANTS																	
LT-Variants	60 - 12	500 - 125	80 - 16	125 - 25	270 - 27	60 - 30	120 - 30	400 - 40	100 - 50	270 - 54	120 - 60	400 - 80	100 - 67	270 - 68	400 - 100	270 - 135	300 - 150
Code No. (REF)	4540002	4540003	4540004	4540006	4540008	4540010	4540012	4540014	4540016	4540018	4540020	4540022	4540024	4540026	4540028	4540032	4540034
Nominal Flow Rate (ml/h)	5	40	5	5	10	2	4	10	2	5	2	5	1.5	4	4	2	2
Nominal Fill Volume (ml)	60	500	80	125	270	60	120	400	100	270	120	400	100	270	400	270	300
Minimal Fill Volume (ml)	30	240	50	60	120	30	60	240	50	135	60	240	50	135	240	135	150
Maximal Fill Volume (ml)	65	560	125	125	335	65	125	560	125	335	125	560	125	335	560	335	340
Maximum Retained Volume (ml)	≤2	≤10	≤3	≤3	≤8	≤2	≤3	≤10	≤3	≤8	≤3	≤10	≤3	≤8	≤10	≤8	≤9

STORAGE

TIME TO REACH ROOM TEMPERATURE (APPROX. +23°C ±2°C)																	
Hours From Refrigerator (approx. +2°C to +8°C)	6	12	6	12	12	6	12	12	6	12	12	12	6	12	12	12	12
Hours From Freezer (approx. -18°C)	12	18	12	18	18	12	18	18	12	18	18	18	12	18	18	18	18

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Note: In order to deliver the Nominal Fill Volume, please add the Maximum Retained Volume:
 Example: LT 60-12: 60 ml + 2 ml = 62 ml (Nominal Fill Volume + Maximum Retained Volume = Total Filling Volume)

en Instructions for Use

Description

Easypump® II ST/LT is intended to infuse medication for either continuous, intravenous, subcutaneous or epidural infusion. Chemotherapy, antibiotics, anesthesia and IV pain management are the most common therapies.

For detailed information about the range of drugs typically administered via elastomeric pumps please refer to the Drug Stability List. Medications and/or fluids must be administered per instructions provided by the drug manufacturer.

Easypump® II is a non-electrically driven portable infusion device, enabling the patient to be treated in an ambulatory manner.

Medication is delivered to the patient by positive pressure applied by the elastomeric membrane. The flow rate is determined by the combination of the flow regulation device (flow restrictor) and the positive pressure of the elastomeric membrane. This pressure delivers the solution against the back-pressure of catheters and blood pressure in the veins. Back-pressure affects the flow rate.

Intended purpose

Single use only.

Indication

Easypump® II is designed to give clinicians and nurses the option of delivering pre-determined amounts of medication to the patient in a continuous and accurate manner.

Patient population

Easypump® II is for patients who receive a treatment in an ambulatory setting and/or at home. In general, the device is not restricted to a special gender. Easypump® II can be used for both adult and pediatric patients.

Intended User

Elastomeric pumps should only be applied by healthcare workers who have been adequately trained in this technique.

Contraindications

Elastomeric pumps are contraindicated for the delivery of blood, blood products, insulin, total parenteral nutrition, lipid emulsions and for intra-articular use as well as for life saving drugs (e.g. catecholamines).


Precautions / Cautions

- Temperature dependency: Easypump® is designed to work at room temperature 23 °C +/- 2 °C (73 °F +/- 3.6 °F). The flow restrictor is calibrated to work at 31 °C (88 °F). To maintain a stable flow rate the flow restrictor should be in close contact with the patient's skin at all times (31 °C). For every 1 °C above or below this temperature, the flow rate will increase or decrease by approximately 3%. An increase in temperature results in an increase in flow rate and vice versa.
- If Easypump® II needs to be stored in the refrigerator or freezer, allow the unit to warm up to room temperature before usage. Storage should not exceed a maximum duration of 30 days, deducting the pump specific nominal infusion time (Example: Maximum 30 days – nominal infusion time 5 days = 25 days maximum storage).
- If storage in freezer is intended, line must not be primed since the filter can be damaged.
- MRI compatible.

- During epidural application the flow rate may increase.
- Delivery time can increase significantly as a result of extended storage time.
- Underfilling and /overfilling: Filling the pump less than nominal volume generally results in a shorter delivery time. Filling the pump more than the nominal volume results in a longer delivery time.
- Diluent dependency: Easypump® flow rates are calculated on basis of using 0.9% NaCl. Using dextrose (D5W) as diluent or the addition of any drug of a higher viscosity than normal saline will increase delivery time (e.g. by 10% in case of dextrose [D5W], see figure 1).
- The system should not be primed when the pump is filled with Dextrose, as the flow rate could decrease significantly.
- An uneven inflation of the pump during the filling process has no effect on the flow rate.
- Air bubbles in the silicone sleeve would not affect the functionality of the pump. There is an air vent in the filter which is capable to filter out the entrapped air bubbles. The silicone sleeve has the properties of gas permeability, therefore the air entrapment in the silicone sleeve can be diffused out of the pump.

External pressure such as squeezing or laying on the pump increases the flow rate.

When filled to nominal volume and under the environmental conditions described here, the flow accuracy is +/- 15% of the labeled flow rate for Easypump®.

Warning

- Do not use if package has been opened or is damaged.
- Single use only. Do not re-sterilize or re-use.
- Re-use of single-use devices creates a potential risk for patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness, or death of the patient.
- Elastomeric pumps should not be used in patients with known hypersensitivity to any of the materials of the device.
- Elastomeric pumps should only be applied by healthcare workers who have been adequately trained in this technique.
- There is no alarm or alert when flow interruption occurs which would advise the user of the potential for serious injury or death due to stoppage or under delivery.

Operating Instructions

Directions for filling

Elastomeric pumps are to be filled with reasonable care and only by qualified personnel.

Prior to usage, please inspect for any damage. Don't use if damaged.

Elastomeric pumps can be filled with a Luer Lock syringe or other similar filling devices using an aseptic technique. The tubing should be primed with diluent before adding medication.

Priming (use aseptic technique):

1. Unfold the outer layer properly prior to start the filling process.
2. Close the clamp.
3. Flip open the comfort cap.
4. Twist off the closing cone from the filling port and place it on a sterile surface.
5. Attach the filling device to the filling port.
6. Inject diluent. Repeat as often as necessary.
7. Remove the filling device from the filling port, reattach the closing cone on the filling port and close all caps.
8. Open the closing cone of the patient connector.
9. Open the clamp and prime the system.
10. Close the clamp and reattach the closing cone of the patient connector. Do not re-open the clamp until the pump is connected to the patient.

Specific priming instructions for drugs prone to precipitation (e.g. 5-Fluorouracil):

- Follow steps 1–5 mentioned under "Priming".
- Initially fill Easyump® II with approximately 10 ml of diluent.
- Follow steps 7–10 listed under "Priming".

Note: If not primed correctly, blockage or slow flow can occur due to crystallization in the area highlighted below.

Filling – Adding medication (use aseptic technique):

Always refer to drug manufacturer's package insert for drug reconstitution/dilution, storage procedures and handling.

Repeat steps 2–6 until required volume is achieved.

After adding the medication, flush the filling port with 0.9% sodium chloride solution to avoid high concentration of the drug inside the filling area.

Note: Please make sure, closing cones of filling port and patient connector are properly closed after filling and prior to use!

At completion, the tubing will be filled exclusively with diluent, protecting from precipitation, while the pump reservoir will contain the medication until the infusion is started.

Patient Connection:

1. Open the closing cone of the patient connector.
2. Connect the patient connector of the pump to the patient's access device.
3. Make sure, that the flow restrictor is completely taped on the patient's skin.
4. Make sure, that the filter is not covered by any dressing.
5. Open the clamp for starting the infusion.

Duration of use

See over- and underfilling table in the Appendix

Disposal

The product Easyump II ST/LT is made of raw materials that both in combination and for themselves do not need to be disposed separately. If used correctly according to the Instructions For Use, no thermal or chemical influences can alter the used materials in a manner that would require separate disposal. Disposal of used devices has to be ensured by the healthcare provider to be in accordance to clinical waste management regulation in each respective country where the medical devices are used.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Material

Materials of the elastomeric pumps that constitute the fluid path are: MVQ, PVC (DEHP-free), PMMA, PES, ABS, glass, e-PTFE

Classification

Easyump® elastomeric pumps are manufactured according to ISO 13485 and are in conformity with the following standards:

- European Community (CE Mark) class IIb Medical Device Directive (MDD)
- Product certified by TÜV Product Service No. 0123 in accordance with 93/42/EEC
- ISO 28620 Medical devices – Non-electrically driven portable infusion devices

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B. Braun Melsungen AG