MonoPlus®

BG	Инструкции за употреба	– 2	LT	Naudojimo instrukcija	– 26
CS	Návod k použití	– 4	LV	Lietošanas instrukcija	– 28
DE	Gebrauchsanweisung	– 6	NL	Gebruiksaanwijzing	– 30
DA	Brugsanvisning	– 8	NO	Bruksanvisning	– 32
ET	Kasutusjuhend	– 10	PL	Instrukcja użytkowania	– 34
ES	Instrucciones de uso	– 12	PT	Instruções de utilização	– 36
FI	Käyttöohjeet	– 14	RO	Instrucțiuni de folosire	– 38
FR	Notice	– 16	RU	Инструкция по применению	– 40
EN	Instructions for use	– 18	SV	Bruksanvisning	– 42
EL	Οδηγίες χρήσης	– 20	SL	Navodila za uporabo	– 44
HU	Használati utasítás	– 22	SK	Návod na použitie	– 46
IT	Istruzioni per l'uso	– 24	TR	Kullanım talimatları	– 48

CE 0123

eIFU v.01 506432 - IFU(506888 02-03/22)

B|BRAUN

SHARING EXPERTISE

B. Braun Surgical, S.A.
Ctra de Terrassa, 121
08191 Rubí, Spain

Phone +34 93 586 62 00
Fax +34 93 588 10 96

www.bbraun.com

MonoPlus®

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MonoPlus®

Description of the device

MonoPlus® is a sterile synthetic absorbable monofilament surgical suture made from the homopolymer poly-p-dioxanone. MonoPlus® is coloured violet with the dyestuff D&C violet No. 2 (Colour Index Number 60725), that is approved by the FDA for monofilament sutures. MonoPlus® fulfils all the requirements of the European and the United States Pharmacopoeia for sterile synthetic absorbable monofilament sutures except for diameter that fulfils the requirements of the European Pharmacopoeia. MonoPlus® is an absorbable suture intended for use in soft tissue approximation (general surgery).

Indications for use

MonoPlus® is indicated in cases where an extended soft tissue wound support of more than 4 weeks is desirable. MonoPlus® presents clinical evidence in general surgery.

The patient population to be used MonoPlus® is adult who do not meet any contraindication, there is not enough evidence in paediatric patients, pregnant or breastfeeding women.

Intended user are healthcare professionals that are qualified and familiar with surgical techniques.

Mode of action

When MonoPlus® suture is used, there is a mild inflammatory reaction, which is typical for an endogenous reaction to a foreign body. MonoPlus® degrades in the body by means of hydrolysis to 2-hydroxyethoxyacetic acid, which is subsequently absorbed and eliminated by the body. After implantation, the hydrolytic process leads to a successive decrease of tensile strength and finally to a complete mass degradation. In vitro and in vivo degradation studies indicate that approximately 65 % to 90 % of its original knot pull tensile strength are retained after 28 days. The mass absorption of MonoPlus® is essentially completed after 180 to 220 days.

Contra-indications

MonoPlus® suture material is contraindicated for tissues that require a permanent wound support, for suturing applying high tension to the tissue or for suturing synthetic implants that require a permanent fixation like in cardio-vascular surgery (vascular grafts or cardiac valves). MonoPlus® suture shall not be used for patients with known intolerances or allergies towards any of the components.

Mode of Application

How to proceed in particular depends on the surgical requirements. Sutures shall be selected depending on the wound size, patient condition, suturing technique and surgeon's experience.

In case that a take-off needle is used, after the operation is completed, please grasp securely the end of the thread and hold the needle holder, pulling strand taut, then the needle is released with a straight tug of the needle holder.

Warnings

MonoPlus® must not be re-sterilized. Open unused or damaged packs must be discarded. Do not use MonoPlus® after expiry date.

Do not reuse the product: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

The user shall be familiar with surgical procedures, when using MonoPlus®. The user must take into consideration that the risk of wound dehiscence may vary depending upon the site of application and the type of suture material used.

Consideration should be taken in the use of MonoPlus® in low vascularised tissue as a delayed absorption may occur.

Usage of MonoPlus® may not be advised in case of patients suffering from diseases or conditions which delay the wound healing process. As MonoPlus® is an absorbable suture material, the surgeon should consider the use of supplemental non-absorbable sutures in the closure of the sites which may undergo expansion, stretching or distension, or which may require additional support.

There is no clinical evidence of MonoPlus® on ophthalmic surgery, otosurgery or neurological tissues, therefore its use must be prevented.

MonoPlus® surgical sutures, with needle removed, contain implantable components which are considered Magnetic Resonance (MR) safe. The implantable poly-p-dioxanone sutures are non-metallic and non-conducting materials and therefore, considered Magnetic Resonance (MR) safe.

Precautions

The user shall be familiar with surgical suturing techniques, before employing MonoPlus® suture materials. When working with suture materials great care shall be taken to ensure that the use of surgical instruments, such as tweezers and needle holders, does not cause the material to be damaged by being pinched or kinked. Care shall be taken to avoid damaging the needle when using the suture material. Always grasp the needle in a section 1/3 to 1/2 of the distance from the fiber attachment end to the needle point, never at the end where the fibre is attached or the needle point. Grasping the needle at the area of its point could impair the penetration performance and cause a fracture of the needle.

Grasping the needle close to the fibre attachment end could cause bending, breakage. Reshaping needles must be avoided and may result in a loss of their strength and resistance towards bending and breaking. Exertion of inappropriate tension on the needle-fiber fixation can lead to a detachment of the fiber from the needle. Care must also be taken to avoid needle stick injuries when handling with surgical needles.

Care shall be taken that the knots are positioned properly, and adequate knot security is given.

In certain applications (e. g. in orthopaedics) immobilization by external support may be employed besides the surgical use of MonoPlus®.

Open the sterile pouch with caution. Avoid contact of the inner sterile part with the outer non-sterile part of the packaging.

MonoPlus® is a single use device. Do not re-use MonoPlus®

Side effects

As with all other suture materials long contact with salt solutions, such as urine and bile, can lead to lithiasis.

As for any other sutures, after implantation the following side effects may occur occasionally: transient local irritation or inflammatory foreign body reaction. An occasional stitch sinus or abscess, infection at the wound site, seroma, haematoma, pain, suture extrusion, suture granuloma and wound dehiscence (causing hernia or burst abdomen) may not be excluded. Existing infections may occasionally be enhanced by any foreign body.

As for all absorbable sutures, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing and/or delayed absorption in tissue with poor blood supply may not be excluded.

Sterilization

MonoPlus® sutures are sterilized by ethylene oxide gas.

The absence of channels along the sealing of the sterile barrier system must be visually checked as well as the absence of perforations on it. In case of identifying such defects, the device shall be discarded in the recommended manner.

Storage

No special storage conditions are required.


Device Disposal

Once the surgical procedure is completed, the different components of MonoPlus® have to be discarded in special intended containers. It is the responsibility of the user to determine if disposal material is hazardous according to federal, state and local regulations. Dispose of contents and container to comply with applicable local, state, national and international regulation. Place all needles in a sharps disposal container immediately after they have been used. Dispose of used sharps disposal containers according to your community guidelines.

Effective

Effective

Symbols used on labelling

	Do not re-use
	Use-by date
	Batch code
	Do not re-sterilize
	Violet
	Polydioxanone
	Dyed / Monofilament / Absorbable
	Sterilized using ethylene oxide
	Single sterile barrier system with protective packaging inside
	Do not use if package is damaged and consult instructions for use
	Date of manufacture
	Consult instructions for use or consult electronic instructions for use
	Medical device
	Catalogue number
	Unique Device Identifier

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UDI

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Title: eIFU MONOPLUS Initiator: Susana ? Calvillo

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Calvillo, Susana (calvsues)
Title: SUT-PR-ES04-Técnico Sistemas Impresión
Date: Wednesday, 06 April 2022, 14:53 W. Europe Daylight Time
Meaning: Document signed as Author

UserName: García, Guillermo (garcgues)
Title: SU-PR-ES04-Technical support manager
Date: Thursday, 07 April 2022, 10:25 W. Europe Daylight Time
Meaning: Precheck of Document

UserName: Orús, Silvia Salut (orussies)
Title: Global Regulatory Affairs Manager CoE CT
Date: Thursday, 07 April 2022, 13:19 W. Europe Daylight Time
Meaning: Approve Document