



B. BRAUN VET CARE

Manuplast[®] VET

EXAMINATION GLOVES FOR
LARGE ANIMALS

Manuplast® VET

Soft, smooth and strong


- High softness due to the ingredient EVA
- Smooth surface – easy to slip on and off
- Strong and proven quality

- Polyethylene (PE), ethylene vinyl acetate (EVA)
- Free of powder, latex proteins and plasticisers
- Flat, ambidextrous shape, welded seam
- Quality tests according to:
ISO 527 (tensile strength, elongation at break)
ISO 7765 (impact resistance)
ISO 6383 (tear resistance)

Manuplast® Vet PE Gloves

- latex-free, powder-free
- made from polyethylene (PE), ethylene vinyl acetate (EVA)
- gloves for large animals, non-sterile
- arm-length: 90 cm
- wall thickness: 22 µm
- transparent orange

SIZE	SALES UNIT	REF
universal size	100	6075100



B. Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen
Tel. +49 (0) 74 61 95-31012 | info@bbraun-vetcare.de | www.bbraun-vetcare.com