## Manuplast<sup>®</sup> Vet Examination gloves for large animals







## Manuplast<sup>®</sup> Vet Soft, smooth and strong

- High softness due to the ingredient EVA
- Smooth surface easy to slip on and off
- Strong and proven quality

- Polyethylene (PE), ethylene vinyl acetate (EVA)
- Free of powder, latex proteins and plasticisers
- Flat, ambidextrous shape, welded seam
- Quality tests according to: ISO 527 (tensile strength, elongation at break) ISO 7765 (impact resistance) ISO 6383 (tear resistance)

## Manuplast<sup>®</sup> Vet

- Light material for better sensitivity
- Arm-length, 90 cm
- Wall thickness 22 μm
- Transparent orange

universal size

size

sales unit REF 6075100 100 pieces



## Manuplast<sup>®</sup> Vet shoulder

- Improved stability and work simplification due to the neck strap
- Better protection and top level hygiene
- Robust material for good tensile strength
- Shoulder length with neck-strap, 118 cm
- Wall thickness 25 μm
- Transparent green

universal size

size



REF

sales unit

50 pieces