

PICK A  
**COLOUR**  
THE RIGHT THREAD  
FOR EVERY SUTURE



## PICK A COLOUR

### VISIT:

[www.bb Braun-vetcare.com/sutures](http://www.bb Braun-vetcare.com/sutures)  
for information such as videos,  
whitepapers and much more.




# THE VETS CHOICE



ABSORBABLE SUTURE MATERIAL

# Monosyn®



## THE STANDARD IN SOFT TISSUE SURGERY

- Wide indication field
- High knot pull tensile strength
- For universal use in soft tissue surgery and intracutaneous sutures
- Wide range of needle-thread combinations

## INDICATIONS:

- Gynaecology, urology, gastro. tract
- Skin, mucosa
- Sub- and intracutaneous sutures

## FEATURES:



Structure:	Monofilament, uncoated
Chemical composition:	Glyconate
Diameter:	USP 7/0 – USP 2
Colour:	undyed, dyed
Absorption time:	~ 14 days
Knot tensile strength retention:	~ 14 days p. i. 50%
Mass absorption:	~ 60 – 90 days



[www.bbraun-vetcare.com/sutures](http://www.bbraun-vetcare.com/sutures)

B. Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen  
Tel. +49 (0) 7461 95-31012 | [info@bbraun-vetcare.de](mailto:info@bbraun-vetcare.de) | [www.bbraun-vetcare.com](http://www.bbraun-vetcare.com)




# THE SOFT ONE



NON-ABSORBABLE SUTURE MATERIAL

# Dafilon®



## THE SOFT AND ECONOMIC SKIN SUTURE

- Soft and smooth
- Smooth tissue passage

## INDICATIONS:

- Skin closure

## FEATURES:




Structure:	Monofilament, uncoated
Chemical composition:	Polyamide (Nylon)
Diameter:	USP 11/0 – USP 5
Colour:	blue, black
Knot tensile strength retention:	non-absorbable



[www.bbraun-vetcare.com/sutures](http://www.bbraun-vetcare.com/sutures)


B. Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen  
Tel. +49 (0) 7461 95-31012 | [info@bbraun-vetcare.de](mailto:info@bbraun-vetcare.de) | [www.bbraun-vetcare.com](http://www.bbraun-vetcare.com)



THE  
RESILIENT ONE

ABSORBABLE SUTURE MATERIAL

# MonoPlus®



## WHEN SUTURING REQUIRES TIME

- Degradation profile designed where an extended wound support is desirable.
- Wide range of needle-thread combinations
- Good knot pull tensile strength retention

## INDICATIONS:

- Orthopaedic surgery, capsular suturing
- Colon surgery
- Abdominal wall closure

## FEATURES:



Structure:	Monofilament, uncoated
Chemical composition:	Poly-p-dioxanone
Diameter:	USP 7/0 – USP 2
Colour:	violet
Absorption time:	~ 35 days
Knot tensile strength retention:	~ 35 days p. i. 50%
Mass absorption:	~ 180 – 210 days



[www.bbraun-vetcare.com/sutures](http://www.bbraun-vetcare.com/sutures)

B. Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen  
Tel. +49 (0) 7461 95-31012 | [info@bbraun-vetcare.de](mailto:info@bbraun-vetcare.de) | [www.bbraun-vetcare.com](http://www.bbraun-vetcare.com)





THE FLEXIBLE ONE

NON-ABSORBABLE SUTURE MATERIAL

# Optilene®



## FLEXIBILITY FOR SPECIAL COMFORT

- Soft and smooth
- Good elasticity and elongation properties
- Proven reliable tensile strength
- Controlled knot run-down

## INDICATIONS:

- Hernias and vessel surgery
- Tendons and capsular suturing
- Skin closure

## FEATURES:



Structure:	Monofilament, uncoated
Chemical composition:	Polypropylene/Polyethylene
Diameter:	USP 10/0 – USP 2
Colour:	blue
Knot tensile strength retention:	non-absorbable



[www.bbraun-vetcare.com/sutures](http://www.bbraun-vetcare.com/sutures)

B. Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen  
Tel. +49 (0) 74 61 95-31012 | [info@bbraun-vetcare.de](mailto:info@bbraun-vetcare.de) | [www.bbraun-vetcare.com](http://www.bbraun-vetcare.com)


# THE RELIABLE ONE

**B|BRAUN**  
SHARING EXPERTISE



ABSORBABLE SUTURE MATERIAL

# Novosyn®



## EXCELLENT KNOT FIT FOR RELIABLE SUTURES

- Improved knot pull tensile strength
- Good flexibility performance

## INDICATIONS:

- Ligatures
- Ophthalmology
- Gynaecology, urology, abdominal wall closure

## FEATURES:




Structure:	Polyfilament braided, coated
Chemical composition:	Poly(glycolide-co-l-lactide 90/10)
Coating:	Poly(glycolide-co-l-lactide 35/65) + CaSt
Diameter:	USP 8/0 – USP 2
Colour:	violet
Absorption time:	~ 21 days
Knot tensile strength retention:	~ 21 days p. i. 50%
Mass absorption:	~ 56 – 70 days



[www.bbraun-vetcare.com/sutures](http://www.bbraun-vetcare.com/sutures)

B. Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen  
Tel. +49 (0) 7461 95-31012 | [info@bbraun-vetcare.de](mailto:info@bbraun-vetcare.de) | [www.bbraun-vetcare.com](http://www.bbraun-vetcare.com)



# THE FAST ONE





ABSORBABLE SUTURE MATERIAL

# Monosyn® Quick



## FAST AND MONOFILAMENT

- Easy handling
- Quick absorption
- High knot pull tensile strength
- Good flexibility and elongation

## INDICATIONS:

- Intracutaneous skin closure
- Gynaecology, urology
- Dental surgery

## FEATURES:



Structure:	Monofilament, uncoated
Chemical composition:	Glyconate
Diameter:	USP 6/0 – USP 1
Colour:	undyed
Absorption time:	~ 6 – 7 days
Knot tensile strength retention:	~ 6 – 7 days p. i. 50%
Mass absorption:	~ 56 days



[www.bbraun-vetcare.com/sutures](http://www.bbraun-vetcare.com/sutures)

B. Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen  
Tel. +49 (0) 74 61 95-31012 | [info@bbraun-vetcare.de](mailto:info@bbraun-vetcare.de) | [www.bbraun-vetcare.com](http://www.bbraun-vetcare.com)