









B. BRAUN VET CARE

# THE PERFECT THREAD FOR EVERY SUTURE

OUR B. BRAUN RECOMMENDATIONS







# THE PERFECT THREAD FOR EVERY SUTURE

## OUR B. BRAUN RECOMMENDATIONS

INDICATION	THREAD STRUCTURE	ARTICLE NUMBER	NAME	SALES UNIT
 <b>SOFT TISSUE SURGERY, PARENCHYMATOUS TISSUE</b>	<b>Monofilament, mid-term absorbable</b>	C2022014	Monosyn® violet USP 4/0 (metric: 1.5), length: 70cm, needle: HR22, Racepack	Pack of 36
		G2022014		Pack of 12
		C2022025	Monosyn® violet USP 3/0 (metric: 2), length: 70cm, needle: HR26, Racepack	Pack of 36
		G2022025		Pack of 12
		C2022036	Monosyn® violet USP 2/0 (metric: 3), length: 70cm, needle: HR30, Racepack	Pack of 36
		G2022036		Pack of 12
 <b>INTRACUTANEOUS SUTURE</b>	<b>Monofilament, mid-term absorbable</b>	C2022204	Monosyn® violet USP 4/0 (metric: 1.5), length: 70cm, needle: DS19, Racepack	Pack of 36
		G2022204		Pack of 12
		C2022215	Monosyn® violet USP 3/0 (metric: 2), length: 70cm, needle: DS24, Racepack	Pack of 36
		G2022215		Pack of 12
 <b>LIGATURES</b>	<b>Monofilament, mid-term absorbable</b>	C0058409	Novosyn® violet USP 1 (metric: 4), length: 2 x 70cm	Pack of 36
		G0058409		Pack of 12
 <b>ABDOMINAL WALL / LINEA ALBA</b>	<b>Monofilament, long-term absorbable</b>	C2024037	MonoPlus® violet USP 0 (metric: 3.5), length: 70cm, needle: HR30, Racepack	Pack of 36
		G2024037		Pack of 12
		C2024058	MonoPlus® violet USP 1 (metric: 4), length: 70cm, needle: HR37S, Racepack	Pack of 36
		G2024058		Pack of 12
 <b>SKIN NON-ABSORBABLE</b>	<b>Monofilament, non-absorbable</b>	C0935204	Dafilon® blue USP 4/0 (metric: 1.5), length: 75cm, needle: DS19	Pack of 36
		G0935204		Pack of 12
		C0935352	Dafilon® blue USP 3/0 (metric: 2), length: 75cm, needle: DS24	Pack of 36
		G0935352		Pack of 12
 <b>CALLUS</b>	<b>Monofilament, non-absorbable</b>	G0068615	Novosyn® violet USP 8/0 (metric: 0.4), length: 15cm, needle: DLM6	Pack of 12

# THE PERFECT THREAD FOR EVERY SUTURE


## OUR B. BRAUN RECOMMENDATIONS

INDICATION	THREAD STRUCTURE	ARTICLE NUMBER	NAME	SALES UNIT
 <b>LASHLINE</b>	<b>Monofilament, mid-term absorbable</b>	C2022204	Monosyn® violet USP 4/0 (metric: 1.5), length: 70cm, needle: DS19	Pack of 36
		G2022204		Pack of 12
 <b>HERNIAS</b>	<b>Monofilament, non- absorbable</b>	C3090041	Optilene® USP 3/0 (metric: 2), length: 75cm, needle: HR26, Racepack	Pack of 36
		G3090041		Pack of 12
 <b>OTHEMATOMA</b>	<b>Monofilament, non- absorbable</b>	C0932655	Dafilon® blue USP 3/0 (metric: 2), length: 75cm, needle: GS60	Pack of 36
		G0932655		Pack of 12
 <b>JOINT CAPSULE</b>	<b>Monofilament, long-term absorbable</b>	C0024131	MonoPlus® violet USP 2/0 (metric: 3), length: 70cm, needle: HR17S	Pack of 36
 <b>DENTAL</b>	<b>Monofilament, short-term absorbable</b>	C0025052	Monosyn® Quick undyed USP 4/0 (metric: 1.5), length: 70cm, needle: HR17	Pack of 36
 <b>SOFT TISSUE SURGERY, PARENCHYMATOUS TISSUE &amp; DENTAL</b>	<b>Monofilament, mid-term absorbable</b>	C2022004	Monosyn® violet USP 4/0 (metric: 1.5), length: 70cm, needle: HR17, Racepack	Pack of 36
		G2022004		Pack of 12

# THE PERFECT THREAD FOR EVERY SUTURE


## OUR B. BRAUN RECOMMENDATIONS

### NEEDLE CURVATURES



- S 1/8 circle
- V 1/4 circle
- D 3/8 circle
- H 1/2 circle
- F 5/8 circle

### MORE NEEDLE CURVATURES



### DESCRIPTION

Needle curvature	Needle shape	Needle tip	Needle length	Needle size
S 1/8 circle	G Cutting, golf ball technology	T Trocar tip	Distance between reinforcement zone and needle tip, dimensions of long end in mm	f Thin-bodied needle
V 1/4 circle	R Round body	N Blunt tip		ff Extra-thin-bodied needle
D 3/8 circle	S Cutting	S Cutting tip (inside)		s Strong needle body
H 1/2 circle	L Lancet	C Short cutting tip		ss Very strong Needle body
F 5/8 circle		M Microtip		v Break-off needle
G Straight		MP Micro Precision Tip		b Black needle
P Asymptotic		m Micro Needle		TO Take-Off removable needle. With a straight tug of the needle holder, the needle can be removed from the suture fiber
J Hook		mV Microvascular Needle		
SK Ski needle				

### Example

HRT 37s = 1/2 circle round-bodied needle with trocar tip, length 37 mm, strong needle body

**CONTACT** | B.Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen

Tel. +49 (0) 7461 95-2191 | Fax +49 (0) 7461 71512 | [bestellung\\_vetcare@aesculap.de](mailto:bestellung_vetcare@aesculap.de) | [www.bbraun-vetcare.com](http://www.bbraun-vetcare.com)