

Askina® Skin Freeze Cryotherapeutic device for treating benign skin lesions

TECHNICAL DATA SHEET



Device variants		Askina® Skin Freeze (2 mm applicator)	Askina® Skin Freeze (5 mm applicator)
Kit components		150 mL aerosol can, box with 50 applicators, instructions for use	
REF		9380701	9380809
Applicator sizes	tip	2 mm tip for standard therapy	5 mm tip for small lesions
	length	150 mm	150 mm
Logistics	packaging	kit package for 50 treatments, carton with 10 kits	
	shelf life	2 years	
Quality certificate		ISO 9001	
Product classification		Medical device acc. to MDD 93/42/EEC, CE class II a	

Askina® Skin Freeze Cryotherapeutic device for treating benign skin lesions

TECHNICAL DATA SHEET

Material Data

150 mL aerosol can



can	aluminum
spraying device	polyethylene
protective cap	polypropylene
applicator carrier	polypropylene


Cryogenic liquid



For hazard and safety information:
please request safety data sheet.

dimethyl ether	(95%)
isobutane	(3%)
propane	(2%)

Applicator



stick	polypropylene
tip	cotton bead, polyurethane foam cover