

Askina® Skin Freeze

CRYOTHERAPEUTIC KIT FOR PROFESSIONAL USE | TECHNICAL DATA



PRODUCT CLASSIFICATION

Active medical device, CE class IIa

QUALITY CERTIFICATES

EN ISO 13485:2016 (CE 0344)

MANUFACTURER



Koninklijke Utermöhlen NV
De Overweg 1, 8471 ZA Wolvega, The Netherlands

DISTRIBUTOR



B. Braun Melsungen AG
Carl-Braun-Str. 1, 34212 Melsungen
Germany

Askina[®] Skin Freeze

CRYOTHERAPEUTIC KIT FOR PROFESSIONAL USE | TECHNICAL DATA

PORTABLE CRYOSURGICAL SYSTEM FOR THE REMOVAL OF WARTS AND OTHER SKIN LESIONS

KIT COMPONENTS

aerosol can filled with 170 ml cryogenic liquid
box with applicators
instructions for use

DEVICE VARIANTS

REF	Applicator sizes	Applicators/pack
9380702	2 mm tip for small lesions	60
9380703	5 mm tip for standard therapy	50
9380704	Mix of 2 and 5 mm applicators	2 x 30

LOGISTICS

packaging kit package for 50–60 treatments, carton with 6 kits
shelf life 3 years
storage conditions



MATERIAL DATA




170 ml aerosol can
can aluminum
spraying device polyethylene
protective cap polypropylene



For hazard and safety information: please request safety data sheet.

Cryogenic liquid DANGER

dimethyl-ether, propane, butane and iso-butane
Extremely flammable aerosol. Pressurized container: May burst if heated. Keep out of reach of children. Do not breathe spray. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not pierce or burn, even after use. Do not spray on an open flame or other ignition source. Protect from sunlight. Do not expose to temperatures exceeding 50°C. Dispose of contents/container in accordance with local regulations.



Applicators stick tip

free of natural rubber latex
polypropylene
cotton bead, polyurethane foam cover