Aesculap Scope Container

Think Sterile Supply. Think Aesculap.



Aesculap Sterile Technology



Suitable sterilisation processes

- Steam sterilisation: fractionated vacuum process
- Low-temperature process: hydrogen peroxide

plasma sterilisation in STERRAD[®] 100S, STERRAD[®] 200, STERRAD[®] NX, STERRAD[®] 100NX

Not for gravitational process!

Separate storage and provision of the scopes

 Separate, quick and reliable aseptic supply, independent of instrument sets

STERRAD* is a registered trademark of Advanced Sterilisation Products (Johnson & Johnson company)

 Protection of sensitive scopes including optical light cable during reprocessing, transportation and storage.



Individual design with the Aesculap organisation system



The new gasket can easily be replaced when necessary






Protected perforation field



The handle on the tray lid allows simple aseptic removal




stackable
tamper proof seal





Lid JM020	LxWxH mm 539 x 152 x 28	for the hydrogen peroxide plasma sterilisation in STERRAD® 100S, STERRAD® 200, STERRAD® NX, STERRAD® 100NX
Body JM021	LxWxH mm 540 x 144 x 65	for the hydrogen peroxide plasma sterilisation in STERRAD® 100S, STERRAD® 200, STERRAD® NX, STERRAD® 100NX
JF167	LxWxH mm 235 x 118	Disposable filter for the hydrogen peroxide plasma sterilisation in STERRAD® 100S, STERRAD® 200, STERRAD® NX, STERRAD® 100NX Order unit (PAK): Pack of 100
JK098	LxWxH mm 245 x 100 x 6	Universal filter holder







			Markin	g pen for Aesculap reusable filters	
JK086				nit (PAK): Pack of 3	
	Colour	Outsid	e dim.	Identification label with	
JG744	white	80 x 35 80 x 35 80 x 35		adhesive back and process indicator for steam sterilisation.	
JG745 JG746	red blue			Order unit (PAK): Pack of 1000	
JG747	green	80 x 3	5		
JG748	yellow	80 x 3	5		
	Outside	dim.		cation label without adhesive	
JG743	80 x 35		back and process indicator for ste sterilisation.		
			Order u	nit (PAK): Pack of 1000	
	Outside			ification label with adhesive back	
10	0705		and pro		
JG750	67 x 35		plasma	cess indicator for STERRAD [®] gas sterilisation.	
JG750	67 X 35			-	
JG750	Outside	dim.	Order u	sterilisation.	
JG750 JG751		dim.	Order u Identifi back ar	sterilisation. nit (PAK): Pack of 100 cation label with adhesive nd process indicator for steam	
	Outside	dim.	Order u Identifi back ar sterilisa	sterilisation. nit (PAK): Pack of 100 cation label with adhesive	

Seals

JG738 JG739

JG740



Seal - safeguarding the container against unauthorized opening.
Order unit (PAK): Pack of 100
Order unit (PAK): Pack of 1000

Repair stamp – for marking instruments for repair or containers holding materials due for repair. *Order unit (PAK): Pack of 1000*



Identification	labola
Identification	laueis

JG646B	blue	Identification labels with lettering one line
JG647B	green	or two lines (character height 5 mm), up to
JG648B	gold	13 characters per line
JG650B	black	Dimensions: 50 x 18 mm
JG652B	olive-brown	Note: Process aluminum identification
JG653B	violet	labels with pH-neutral cleaning agent.
JG653B		labels with pH-neutral cleaning agent.





JG645U	red
JG646U	blue
JG647U	green
JG648U	gold
JG650U	black
JG652U	olive-brown
JG653U	violet
JG654U	orange

Specifications as above

Note: Process aluminum identification labels with pH-neutral cleaning agent.

Outside dim. JG736 56 x 18	Bar code labels, plastic, 128A font, max. code length 10 characters, max. lettering length 10 characters. Special note: Please specify the required designation in your order.
Outside dim. IC-V006 56 x 18	Bar code labels, plastic, for instacount [®] , 128A font, max. code length 10 characters, max. lettering length 20 characters. Special note: Please specify the required designation and the serial number of the instacount set in your order.
Organisation System	Eyepiece holder
JG367	-
JG368	Sleeve holder
JG300	Fixation pin Order unit (PAK): Pack of 10
JG302	Screwdriver





Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany Phone +49 074 61 95-0 | Fax +49 074 61 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

All rights reserved. Technical alterations are possible. This leaflet may be used in no way whatsoever other than for the offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the right to recall our catalogues and pricelists and to take legal action.

Brochure No. C26502