



B. BRAUN VET CARE

AESCLAP INSTRUMENTS

SET CONFIGURATION AND STANDARD INSTRUMENTS FOR VETERINARY MEDICINE

BASIC SURGERY SET | VF231R

Art-No.	Indication	Piece
BB073R	Scalpel handle #3, 125 mm	1
BB084R	Scalpel handle #4, 135 mm	1
BD557R	Standard Tissue Forceps, 1 x 2, 145 mm	1
BC324R	Surgical Scissors, straight, s/b, 145 mm	1
BC605R	Metzenbaum Scissors, curved, 145 mm	1
BD700R	Adson-Brown Tissue Forceps, 120 mm	1
BH104R	Hartmann Mosquito forceps, delicate, straight, 100 mm	2
BH120R	Halsted-Mosquito forceps, delicate, straight, 1x2, 125 mm	2
BH425R	Pean Forceps, curved, 140 mm	2
BM065R	TC Mayo-Hegar Needle Holder, 150 mm	1
FK630R	Volkman Bone Curette, 8 mm, 170 mm	1



SOFT TISSUE SURGERY SET | VF230R

Art.-No.	Indication	Piece
BB073R	Scalpel handle #3, 125 mm	1
BB084R	Scalpel handle #4, 135 mm	1
BC011R	Stevens scissors, curved, blunt/blunt, 115 mm	1
BC324R	Surgical scissors, straight, s/b, 145 mm	1
BC605R	Metzenbaum scissors, curved, 145 mm	1
BD511R	Adson Tissue Forceps, fine, 1 x 2 t, 120 mm	1
BD557R	Surgical scissors, 1 x 2 t, 145 mm	1
BD700R	Adson-Brown Tissue Forceps, 120 mm	1
BF431R	Backhaus Towel Clamp, 90 mm	5
BH104R	Hartmann Mosquito forceps, delicate, straight, 100 mm	2
BH120R	Halsted-Mosquito forceps, delicate, straight, 1 x 2, 125 mm	2
BH425R	Pean forceps, curved, 140 mm	2
BM065R	TC Mayo-Hegar needle holder, 150 mm	1
BM806R	Deschamps Needle, blunt, 215 mm	1
BT008R	Senn-Miller Rector., sh., 8 x 7/18 x 5.5 mm	1
BV074R	Weitlaner retractor, 3 x 4 t., blunt, 130 mm	1
EA016R	Allis forceps, 5 x 6, 155 mm	1
FK321R	Elevators-Holder, curved., round-tip 6 mm, 190 mm	1
FK630R	Volkman Bone Curette, 8 mm, 170 mm	1
VF225R	Castration hook for cats, 170 mm	1



CASTRATION SET FOR CATS | VF232R

Art.-No.	Indication	Piece
BB073R	Scalpel handle #3, 125 mm	1
BC605R	Metzenbaum scissors, curved, 145 mm	1
BD700R	Adson-Brown Tissue Forceps, 120 mm	1
BH120R	Halsted-Mosquito forceps, delicate, straight, 1 x 2, 125 mm	2
BM065R	TC Mayo-Hegar needle holder, 150 mm	1
VF225R	Castration hook f. cats, 170 mm	1



STORAGE SYSTEM FOR INSTRUMENTS AND SETS

Abbildung

Art.-No. Indication

STORAGE SYSTEM FOR BASIC SURGERY SET




JN092 Veterinary Sterile Container, 312 x 190 x 46 mm



JG381R Instrument tray w. lid, 273 x 176 x 30 mm

STORAGE SYSTEM FOR SOFT TISSUE SURGERY SET



JG386R Instrument tray (empty)



JN094 Veterinary Sterile Container, 313 x 190 x 65 mm


Abbildung

Art.-No. Indication

STORAGE SYSTEM FOR CASTRATION SET FOR CATS



JK187 Mini container base, 310 x 140 x 40 mm



JK170 Lid for mini container red




JG387R Instrument tray, 188 x 84 x 30 mm

ORGANIZING ELEMENTS

















JF346 Fixation pin for ring handle instruments








JF945 Cut-to-size silicone mat







SCISSORS

Instrument	Art.-No.	Indication
	BC323R	Surgical scissors, straight, sharp/blunt, 130 mm
	BC423R	Surgical scissors, curved, sharp/blunt, 130 mm
	BC314R	Surgical scissors, straight, blunt/blunt, 145 mm
	BC414R	Surgical scissors, curved, blunt/blunt, 145 mm
	BC260R	TC Metzenbaum scissors, straight, 145 mm
	BC261R	TC Metzenbaum scissors, curved, 145 mm





Instrument	Art.-No.	Indication
	BC601R	Metzenbaum scissors, straight, 145 mm
	BC602R	Metzenbaum scissors, straight, 180 mm
	BC605R	Metzenbaum scissors, curved, 145 mm
	BC606R	Metzenbaum scissors, curved, 180 mm
	BC545R	Mayo scissors, straight, 155 mm
	BC555R	Mayo scissors, curved, 155 mm
	BC166R	Kilner scissors, curved, with flat tips, 125 mm
	BC010R	Stevens scissors, straight, blunt/blunt, 115 mm
	BC011R	Stevens scissors, curved, blunt/blunt, 115 mm







SCISSORS, SCALPEL HANDLE, RETRACTOR

Instrument	Art.-No.	Indication
	BC164R	Strabismus scissors, straight, blunt/blunt, 115 mm
	BC165R	Strabismus scissors, curved, blunt/blunt, 115 mm
	BC294W	TC Suture Wavecut scissors, curved, 145 mm
	BC861R	Lister Bandage scissors, 140 mm
	BC881R	Universal scissors black, 180 mm







Instrument	Art.-No.	Indication
	BB084R	Scalpel handle #4, 135 mm
	BB073R	Scalpel handle #3, 125 mm
	BT008R	Senn-Miller Rector. Sh., 8 x 7/18 x 5.5 mm
	BT272R	Kocher retractor, blunt, 15.5 x 9 mm, 220 mm
	BT359R	Kocher-Langenbeck retractor, 41 x 11 mm, 215 mm
	BT177R	Langenbeck Green retractor, 16 x 6 mm, 160 mm








RETRACTOR, GELPI, FORCEPS

Instrument	Art.-No.	Indication
 A metal surgical retractor with two long, curved arms. Each arm has three small, rounded teeth at the tip. The handles are circular and positioned on the right side.	BV074R	Weitlaner retractor, 3 x 4 t, blunt, 130 mm
 A metal surgical forceps with two curved, pointed jaws. The handles are circular and positioned on the right side.	VB142R	Gelpi small, flat, 93 mm
 A metal surgical forceps with two curved, pointed jaws. The handles are circular and positioned on the right side.	VB143R	Gelpi small, curved, 88 mm
 A metal surgical forceps with two long, curved, pointed jaws. The handles are circular and positioned on the right side.	BV996R	Gelpi Vaginal, 135 mm

Instrument	Art.-No.	Indication
	BH110R	Halsted-Mosquito forceps, delicate, straight, 125 mm
	BH111R	Halsted-Mosquito forceps, delicate, curved, 125 mm
	BH120R	Halsted-Mosquito forceps, delicate, straight, 1 x 2, 125 mm
	BH121R	Halsted-Mosquito forceps, delicate, curved, 1 x 2, 125 mm
	BH104R	Hartmann Mosquito forceps, delicate, straight, 100 mm
	BH105R	Hartmann Mosquito forceps, delicate, curved, 100 mm






FORCEPS, HOOKS, BONE CURETTES

Instrument	Art.-No.	Indication
	BH424R	Pea forceps, straight, 140 mm
	BH425R	Pea forceps, curved, 140 mm
	BH442R	Rochester-Pea forceps, straight, 160 mm
	BH443R	Rochester-Pea forceps, curved, 160 mm
	BH642R	Kocher-Ochsner forceps, straight, 1 x 2, 160 mm
	BH643R	Kocher-Ochsner forceps, curved, 1 x 2, 160 mm

Instrument	Art.-No.	Indication
	EA016R	Allis forceps, 5 x 6, 155 mm
	EA125R	Hartmann intestinal forceps, curved, 200 mm
	VF225R	Castration hook f. cats, 170 mm
	BM806R	Deschamps needle, blunt, 215 mm
	BN104R	Double-Ended Probe, 1.0 mm, 145 mm
	FK842R	Volkman bone curette, double-end, 170 mm
	FK636R	Volkman bone curette #3, 10.0 mm, 170 mm

RASPATORY, ELEVATOR, NEEDLE HOLDER

Instrument	Art.-No.	Indication
	FK300R	Williger bone raspatory, 6.0 mm, 160 mm
	FK321R	Elevators-Handle, curved, round-tip, 6 mm, 190 mm
	FK317R	Elevator Peek-Handle, curved, 3 mm, 185 mm
	FK117R	Bone Elevator, 160 mm, 6 mm
	FK170R	Bone elevator, 10 mm, Wi. 220 mm


Instrument	Art.-No.	Indication
	BM065R	TC Mayo-Hegar needle holder, 150 mm
	BM235R	Mayo-Hegar Needle Holder, 150 mm
	BM151R	TC Mathieu Needle Holder, 170 mm
	BM361R	Mathieu Needle Holder, 170 mm
	BM128R	TC Hegar-Olsen Needle Holder/Scissors, 145 mm

AESCULAP INSTRUMENTS


SUPPORTING YOUR PROCESS:



V760 Container




V066 Sutures



V087 Osteosynthesis





AESCULAP® – a B. Braun brand

B. Braun Vet Care GmbH | Am Aesculap-Platz | 78532 Tuttlingen
Fax +49 (0) 74 61 715 12 | info@bbraun-vetcare.com | www.bbraun-vetcare.com

The main product trademark 'Aesculap' and the product trademarks 'Durotip' and 'Durogrip' are registered trademarks of Aesculap AG. Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.