

Caresite – Needlefree Luer Access Device

For all injections and infusions



New!

Placeholder
specialised dealer

B | BRAUN
SHARING EXPERTISE

Caresite – Needlefree Luer Access Device

For all injections and infusions

Special offer



Easy to visually inspect



Easy to clean



Easy to connect

Your advantages:

Helps to prevent catheter occlusion





- Positive displacement valve design

Safe and hygienic conditions

- needlefree injection and aspiration
- smooth surface supports easy disinfection
- visible fluid path

Convenient and easy to use

- Therapy interruptions at anytime without additional cones etc.

	Description	Art.-No.	Length in cm	Priming Volume	Inner-Ø in mm	Min. quantity
	Caresite	415122-01	-	0.22 ml	-	10
	Caresite with extension and Spin-Lock connector	470108-01	15	0.9 ml	2.8	10
	2 Caresite with Y-extension and Spin-Lock	470106-01	18	0.9 ml	1.3	10
	Intrafix Primeline Y-Caresite	4062158C	180	-	3	10

Discount

Hospital

Street

Zip-Code, City

Date, Signature

Placeholder specialised dealer

B | BRAUN
SHARING EXPERTISE

B. Braun Vet Care GmbH
info@bbraun-vetcare.de
www.bbraun-vetcare.com