### Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 Issue date: 1/31/2023 Revision date: 1/31/2023 Supersedes: 7/8/2020 Version: 2.1 SDS No: 00056-0072

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form Product name UFI

: Mixture : Softasept N coloured / uncoloured 2R0V-C762-300Q-75YR VT0V-U7VF-E006-VHJT

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Use of the substance/mixture

: Alcoholic cutaneous solution

## 1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Manufacturer	Supplier
B. Braun Medical AG	B. Braun Melsungen AG
Seesatz 17	Carl-Braun-Straße 1
CH-6204 Sempach	34212 Melsungen – Deutschland
Switzerland	Zentrale Service-Bereich / Logistik und Supply Chain
T +41 (0) 58 / 258 50 00	T +49 (0) 5661 / 71-4422
info.bbmch@bbraun.com	logistics.service@bbraun.com
E-mail address of competent person responsible for the SDS: sds@gbk-	ingelheim.de

#### 1.4. Emergency telephone number

Emergency number

: INTERNATIONAL: +49 - (0) 6132 - 84463, GBK GmbH (24h - 7d/w - 365d/a)

#### **SECTION 2: Hazards identification**

#### 2.1. Classification of the substance or mixture

<b>Classification according to Regulation</b>	(EC) No.	1272/2008 [CLP]
---	----------	-----------------

Flammable liquids, Category 2	H225
Serious eye damage/eye irritation, Category 2	H319
Full text of H- and EUH-statements: see section 16	

#### Adverse physicochemical, human health and environmental effects

Highly flammable liquid and vapour. Causes serious eye irritation.

#### 2.2. Label elements

Labelling according	a to Regulation	(FC) No. 12	72/2008 [CLP]
Labening according	g to negulation		

Laboling according to regulation (L	
Hazard pictograms (CLP)	
Signal word (CLP)	: Danger
Hazard statements (CLP)	: H225 - Hig
	H319 - Cau

Precautionary statements (CLP)




- hly flammable liquid and vapour.
- uses serious eye irritation.
- : P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P260 - Do not breathe vapours.
- P280 Wear eye protection.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

# B BRAUN

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072



#### 2.3. Other hazards

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII Contains no PBT/vPvB substances  $\geq$  0.1% assessed in accordance with REACH Annex XIII

The mixture does not contain substance(s) included in the list established in accordance with Article 59(1) of REACH for having endocrine disrupting properties, or is not identified as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at a concentration equal to or greater than 0,1 %

### **SECTION 3: Composition/information on ingredients**

#### 3.1. Substances

Not applicable

#### 3.2. Mixtures

#### Comments

: Alcoholic solution

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Ethanol	CAS-No.: 64-17-5 EC-No.: 200-578-6 EC Index-No.: 603-002-00-5 REACH-no: 01-2119457610- 43	< 75	Flam. Liq. 2, H225 Eye Irrit. 2, H319
Propan-2-ol	CAS-No.: 67-63-0 EC-No.: 200-661-7 EC Index-No.: 603-117-00-0 REACH-no: 01-2119457558- 25	< 15	Flam. Liq. 2, H225 Eye Irrit. 2, H319 STOT SE 3, H336

Specific concentration limits:		
Name	Product identifier	Specific concentration limits
	CAS-No.: 64-17-5 EC-No.: 200-578-6 EC Index-No.: 603-002-00-5 REACH-no: 01-2119457610- 43	( 50 ≤C ≤ 100) Eye Irrit. 2, H319

Full text of H- and EUH-statements: see section 16

### Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

### **SECTION 4: First aid measures**

#### 4.1. Description of first aid measures

First-aid measures general	: Data of item 4 do partly not refer to the use and the regular employing of the product (in this sense consult package leaflet and expert information), but to liberation of major amounts in case of accidents and irregularities.
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing.
First-aid measures after skin contact	: Rinse skin with water/shower. Take off immediately all contaminated clothing.
First-aid measures after eye contact	: Wash immediately with plenty water (during 20 minutes), also under eyelids. Remove contact lenses, if present and easy to do. Continue rinsing. Call a physician immediately.
First-aid measures after ingestion	: Do not induce vomiting. Call a physician immediately. Drink plenty of water. Never give anything by mouth to an unconscious person.
4.2. Most important symptoms and eff	fects, both acute and delayed
Symptoms/effects after inhalation	: Narcotic in high concentrations.
Symptoms/effects after skin contact	: Repeated or prolonged skin contact may cause dermatitis and defatting.
Symptoms/effects after eye contact	: Eye irritation.
Symptoms/effects after ingestion	: Ingestion causes irritation of upper respiratory system and gastrointestinal disturbance.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures	
5.1. Extinguishing media	
Suitable extinguishing media Unsuitable extinguishing media	<ul><li>Water spray. Dry powder. Carbon dioxide. Alcohol-resistant foam.</li><li>high volume water jet.</li></ul>
5.2. Special hazards arising from the subst	tance or mixture
Fire hazard Explosion hazard Hazardous decomposition products in case of fire	<ul> <li>Highly flammable liquid and vapour.</li> <li>Product is not explosive. Explosive vapour/air mixtures may be formed.</li> <li>Toxic fumes may be released. Carbon oxides (CO, CO2).</li> </ul>
5.3. Advice for firefighters	
Precautionary measures fire Firefighting instructions Protection during firefighting Other information	<ul> <li>Cool endangered containers with water spray jet.</li> <li>Fight fire from safe distance and protected location.</li> <li>Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.</li> <li>Vapours are heavier than air and may spread along floors. Collect contaminated fire-fighting water, avoid any release into the sewerage. Fire residues and contaminated firefighting water must be disposed of in accordance with the local regulations.</li> </ul>

SECTION 6: Accidental release measures		
6.1. Personal precautions, protecti	ve equipment and emergency procedures	
6.1.1. For non-emergency personnel		
Emergency procedures	: Ventilate spillage area. No open flames, no sparks, and no smoking. Avoid contact with skin and eyes.	
6.1.2. For emergency responders		
Protective equipment	: Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".	
6.2. Environmental precautions		

Avoid release to the environment.

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

6.3. Methods and material for containment and cleaning up		
For containment	: Dike and contain spill.	
Methods for cleaning up	: Take up liquid spill into absorbent material. Notify authorities if product enters sewers or public waters.	
Other information	: Dispose of materials or solid residues at an authorized site.	
6.4. Reference to other sections		

Refer to protective measures listed in sections 7 and 8. For further information refer to section 13.

SECTION 7: Handling and storage	
7.1. Precautions for safe handling	
Precautions for safe handling	: Ensure good ventilation of the work station. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Ground/bond container and receiving equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Flammable vapours may accumulate in the container. Use explosion-proof equipment. Wear personal protective equipment. Avoid contact with eyes.
Hygiene measures	: Do not eat, drink or smoke when using this product. Wash contaminated clothing before reuse. Always wash hands after handling the product.
7.2. Conditions for safe storage, includin	g any incompatibilities
Technical measures Storage conditions Incompatible materials Information on mixed storage	<ul> <li>Ground/bond container and receiving equipment.</li> <li>Store in a well-ventilated place. Keep cool. Keep container tightly closed.</li> <li>oxidizing materials. alkali metals. alkaline earth metals.</li> <li>Keep away from food, drink and animal feeding stuffs.</li> </ul>

#### 7.3. Specific end use(s)

See Section 1.

## SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

#### 8.1.1 National occupational exposure and biological limit values

No additional information available

#### 8.1.2. Recommended monitoring procedures

Monitoring methods	
Monitoring methods	A specific exposure sampling method is not available.
Biological monitoring methods	A specific exposure sampling method is not available

#### 8.1.3. Air contaminants formed

No additional information available

#### 8.1.4. DNEL and PNEC

Ethanol (64-17-5)	
DNEL/DMEL (General population)	
Acute - local effects, inhalation 1900 mg/m <sup>3</sup>	
Long-term - systemic effects,oral	343 mg/kg bodyweight/day
Long-term - systemic effects, inhalation	950 mg/m³
Long-term - systemic effects, dermal	343 mg/kg bodyweight/day

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

Ethanol (64-17-5)			
PNEC (Water)			
PNEC aqua (freshwater)	0.96 mg/l		
PNEC aqua (marine water)	0.79 mg/l		
PNEC (Sediment)			
PNEC sediment (freshwater)	3.6 mg/kg dwt		
PNEC sediment (marine water)	2.9 mg/kg dwt		
PNEC (Soil)			
PNEC soil	0.63 mg/kg dwt		
PNEC (STP)			
PNEC sewage treatment plant	580 mg/l		
Propan-2-ol (67-63-0)			
DNEL/DMEL (Workers)			
Long-term - systemic effects, dermal	888 mg/kg bodyweight/day		
Long-term - systemic effects, inhalation	500 mg/m³		
DNEL/DMEL (General population)			
Long-term - systemic effects,oral	26 mg/kg bodyweight/day		
Long-term - systemic effects, inhalation	88 mg/m <sup>3</sup>		
Long-term - systemic effects, dermal	319 mg/kg bodyweight/day		
PNEC (Water)	PNEC (Water)		
PNEC aqua (freshwater)	140.9 mg/l		
PNEC aqua (marine water)	140.9 mg/l		
PNEC (Sediment)			
PNEC sediment (freshwater)	552 mg/kg dwt		
PNEC sediment (marine water)	552 mg/kg dwt		
PNEC (Soil)			
PNEC soil	28 mg/kg dwt		

#### 8.1.5. Control banding

No additional information available

#### 8.2. Exposure controls

#### 8.2.1. Appropriate engineering controls

#### Appropriate engineering controls:

Ensure good ventilation of the work station.

#### 8.2.2. Personal protection equipment

#### Personal protective equipment:

Data of item 8 do partly not refer to the use and the regular employing of the product (in this sense consult package leaflet and expert information), but to liberation of major amounts in case of accidents and irregularities.

#### 8.2.2.1. Eye and face protection

#### Eye protection:

Eyewash bottle with clean water (EN 15154)

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

Eye protection			
Туре	Field of application	Characteristics	Standard
Safety glasses with side shields	Liquid splashes may occur		EN 166

#### 8.2.2.2. Skin protection

#### Hand protection:

This recommendation refers exclusively to the chemical compatibility and the lab test conforming to EN 374 carried out under lab conditions. Requirements can vary as a function of the use. Therefore it is necessary to adhere additionally to the recommendations given by the manufacturer of protective gloves

Hand protection					
Туре	Material	Permeation	Thickness (mm)	Penetration	Standard
Chemically resistant protective gloves	Butyl rubber	6 (> 480 minutes)	0,7		EN ISO 374

#### 8.2.2.3. Respiratory protection

#### **Respiratory protection:**

In case of insufficient ventilation, wear suitable respiratory equipment. No personal breathing protective equipment is normally required

Respiratory protection			
Device	Filter type	Condition	Standard
Breathing apparatus with filter	Type A - High-boiling (>65 °C) organic compounds	In the event of insufficient ventilation:	EN 14387

#### 8.2.2.4. Thermal hazards

No additional information available

#### 8.2.3. Environmental exposure controls

#### Environmental exposure controls:

Avoid release to the environment.

## **SECTION 9: Physical and chemical properties**

#### 9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Colour	: Colourless. red.
Appearance	: Liquid.
Odour	: characteristic.
Odour threshold	: Not available
Melting point	: -114 °C Ethanol
Freezing point	: Not available
Boiling point	: 78 °C Ethanol
Flammability (solid, gas)	: Not applicable
Explosive properties	: Product is not explosive. Explosive vapour/air mixtures may be formed.
Oxidising properties	: Not oxidising.
Explosive limits	: Not available
Lower explosive limit (LEL)	: 3.5 vol % Ethanol
Upper explosive limit (UEL)	: 15 vol % Ethanol
Flash point	: 12 – 14 °C DIN 51755
Auto-ignition temperature	: 425 °C Ethanol
Decomposition temperature	: > 700 °C Ethanol
pH	: 7
Viscosity, kinematic	: Not available
Solubility	: Water: Miscibility
Partition coefficient n-octanol/water (Log Kow)	: Not available

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

Vapour pressure Vapour pressure at 50°C Density	<ul> <li>59 hPa 20°C, Ethanol</li> <li>Not available</li> <li>0.82 - 0.84 g/cm<sup>3</sup> (20°C)</li> </ul>
Relative density	: Not available
Relative vapour density at 20°C	: 2.07 Ethanol
Particle characteristics	: Not applicable

#### 9.2. Other information

#### 9.2.1. Information with regard to physical hazard classes

No additional information available

#### 9.2.2. Other safety characteristics

VOC content Additional information : < 90 % : Solvent content < 90 %

SECTION 10: Stability	v and reactivity
	y and reactivity

#### 10.1. Reactivity

Highly flammable liquid and vapour. No decomposition if stored normally.

#### 10.2. Chemical stability

Stable under normal conditions.

#### 10.3. Possibility of hazardous reactions

Reacts with: oxidizing materials. alkali metals. alkaline earth metals.

#### 10.4. Conditions to avoid

Avoid contact with hot surfaces. Heat. No flames, no sparks. Eliminate all sources of ignition.

#### 10.5. Incompatible materials

Oxidizing agent. alkali metals. alkaline earth metals.

#### 10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced. Thermal decomposition generates : Carbon oxides (CO, CO2).

#### **SECTION 11: Toxicological information**

#### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity (oral) Acute toxicity (dermal) Acute toxicity (inhalation)	<ul> <li>Not classified (Based on available data, the classification criteria are not met)</li> <li>Not classified (Based on available data, the classification criteria are not met)</li> <li>Not classified (Based on available data, the classification criteria are not met)</li> </ul>
Ethanol (64-17-5)	
LD50 oral rat	10470 mg/kg
LD50 oral	6200 mg/kg
LD50 dermal rabbit	15800 mg/kg
LD50 dermal	20000 mg/kg
LC50 Inhalation - Rat (Vapours)	30 mg/l/4h
Propan-2-ol (67-63-0)	
LD50 oral rat	5840 mg/kg

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

Propan-2-ol (67-63-0)	
LD50 dermal rabbit	13900 mg/kg
LD50 dermal	12870 mg/kg
LC50 Inhalation - Rat	> 25 mg/l 4 h
Skin corrosion/irritation	<ul> <li>Not classified (Based on available data, the classification criteria are not met) pH: 7</li> </ul>
Serious eye damage/irritation	: Causes serious eye irritation. pH: 7
Respiratory or skin sensitisation	: Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity	: Not classified (Based on available data, the classification criteria are not met)
Carcinogenicity	: Not classified (Based on available data, the classification criteria are not met)
Reproductive toxicity	: Not classified (Based on available data, the classification criteria are not met)
STOT-single exposure	: Not classified (Based on available data, the classification criteria are not met)
Propan-2-ol (67-63-0)	
STOT-single exposure	May cause drowsiness or dizziness.
STOT-repeated exposure	Not classified (Based on available data, the classification criteria are not met)
Aspiration hazard	: Not classified (Based on available data, the classification criteria are not met)
11.2. Information on other hazards	

#### 11.2.1. Endocrine disrupting properties

No additional information available

#### 11.2.2. Other information

 Potential adverse human health effects and symptoms
 : At high concentrations, the vapours may cause narcosis,Repeated or prolonged exposure may cause skin irritation and dermatitis, due to degreasing properties of the product,Ingestion causes irritation of upper respiratory system and gastrointestinal

## SECTION 12: Ecological information

12.1. Toxicity	
Ecology - general	: The product is not considered harmful to aquatic organisms nor to cause long-term adverse effects in the environment.
Hazardous to the aquatic environment, short-term (acute)	: Not classified (Based on available data, the classification criteria are not met)
Hazardous to the aquatic environment, long-term (chronic)	: Not classified (Based on available data, the classification criteria are not met)
Ethanol (64-17-5)	
LC50 fish 1	11200 mg/l 96 h
EC50 Daphnia 1	5012 mg/l 48 h, Daphnia magna
	275 mg/l

disturbance

ErC50 algae	275 mg/l
Propan-2-ol (67-63-0)	
LC50 fish 1	9640 mg/l Pimephales promelas, 96 h
EC50 Daphnia 1	10000 mg/l Daphnia magna, 48 h
EC50 72h - Algae [1]	1800 mg/l Desmodesmus subspicatus, 72 h

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

12.2. Persistence and degradability		
Ethanol (64-17-5)		
Persistence and degradability	Readily biodegradable.	
Propan-2-ol (67-63-0)		
Persistence and degradability	Readily biodegradable.	
Biodegradation	95 % 21 d, (OECD 301E method)	

### 12.3. Bioaccumulative potential

Ethanol (64-17-5)		
Bioconcentration factor (BCF REACH)	0.06	
Log Pow	-0.3	
Propan-2-ol (67-63-0)		
Log Pow	0.05	

#### 12.4. Mobility in soil

No additional information available

#### 12.5. Results of PBT and vPvB assessment

Softasept N coloured / uncoloured		
This substance/mixture does not meet the PBT cr	iteria of REACH regulation, annex XIII	
This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII		
12.6. Endocrine disrupting properties		
Adverse effects on the environment caused by	• The mixture does not contain substance(s) included in the list established in accordance	

Adverse effects on the environment caused by	: The mixture does not contain substance(s) included in the list established in accordance
endocrine disrupting properties	with Article 59(1) of REACH for having endocrine disrupting properties, or is not identified
	as having endocrine disrupting properties in accordance with the criteria set out in
	Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU)
	2018/605 at a concentration equal to or greater than 0,1 %.

#### 12.7. Other adverse effects

No additional information available

SECTION 13: Disposal consideration	S
13.1. Waste treatment methods	
Waste treatment methods	: Dispose of contents/container in accordance with licensed collector's sorting instructions.
Product/Packaging disposal recommendations	: Empty containers should be taken for recycling, recovery or waste in accordance with local regulation.
Additional information	: Flammable vapours may accumulate in the container.
European List of Waste (LoW) code	: 07 06 04* - other organic solvents, washing liquids and mother liquors

### **SECTION 14: Transport information**

In accordance with ADR / IMDG / IATA / ADN / RID

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

ADR	IMDG	IATA	ADN	RID			
14.1. UN number or ID number							
UN 1987	UN 1987	UN 1987	UN 1987	UN 1987			
14.2. UN proper shippin	I4.2. UN proper shipping name						
ALCOHOLS, N.O.S. (Ethanol ; Propan-2-ol)	ALCOHOLS, N.O.S. (Ethanol ; Propan-2-ol)	Alcohols, n.o.s. (Ethanol ; Propan-2-ol)	ALCOHOLS, N.O.S. (Ethanol ; Propan-2-ol)	ALCOHOLS, N.O.S. (Ethanol ; Propan-2-ol)			
Transport document descr	iption						
UN 1987 ALCOHOLS, N.O.S. (Ethanol ; Propan-2- ol), 3, II, (D/E)	UN 1987 ALCOHOLS, N.O.S. (Ethanol ; Propan-2- ol), 3, II	UN 1987 Alcohols, n.o.s. (Ethanol ; Propan-2-ol), 3, II	UN 1987 ALCOHOLS, N.O.S. (Ethanol ; Propan-2- ol), 3, II	UN 1987 ALCOHOLS, N.O.S. (Ethanol ; Propan-2- ol), 3, II			
14.3. Transport hazard o	class(es)						
3	3	3	3	3			
14.4. Packing group							
II	II	II	II	II			
14.5. Environmental hazards							
Dangerous for the environment: No	Dangerous for the environment: No Marine pollutant: No	Dangerous for the environment: No	Dangerous for the environment: No	Dangerous for the environment: No			
No supplementary information	No supplementary information available						
14.6. Special precaution	s for user						

#### **Overland transport**

Classification code (ADR) Special provisions (ADR) Limited quantities (ADR) Excepted quantities (ADR)	::	F1 274, 601, 640D 1I E2
Packing instructions (ADR)	:	P001, IBC02, R001
Mixed packing provisions (ADR)	:	MP19
Transport category (ADR)	:	2
Hazard identification number (Kemler No.)	:	33
Orange plates	:	33 1987
Tunnel restriction code (ADR)	:	D/E
Transport by sea		
Special provisions (IMDG)	:	274
Limited quantities (IMDG)	:	1 L
Excepted quantities (IMDG)	:	E2
Packing instructions (IMDG)	:	P001
IBC packing instructions (IMDG)	:	IBC02
Tank instructions (IMDG)	:	T7
Tank special provisions (IMDG)	:	TP1, TP8, TP28
EmS-No. (Fire)	:	F-E
EmS-No. (Spillage)	:	S-D
Stowage category (IMDG)	:	В

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

Air transport	
PCA Excepted quantities (IATA)	: E2
PCA Limited quantities (IATA)	: Y341
PCA limited quantity max net quantity (IATA)	: 1L
PCA packing instructions (IATA)	: 353
PCA max net quantity (IATA)	: 5L
CAO packing instructions (IATA)	: 364
CAO max net quantity (IATA)	: 60L
Special provisions (IATA)	: A3, A180
ERG code (IATA)	: 3L
Inland waterway transport	
Classification code (ADN)	: F1
Special provisions (ADN)	: 274, 601, 640D
Limited quantities (ADN)	: 1L
Excepted quantities (ADN)	: E2
Carriage permitted (ADN)	: T
Equipment required (ADN)	: PP, EX, A
Ventilation (ADN)	: VE01
Number of blue cones/lights (ADN)	: 1
Rail transport	
Classification code (RID)	: F1
Special provisions (RID)	: 274, 601, 640D
Limited quantities (RID)	: 1L
Excepted quantities (RID)	: E2
Packing instructions (RID)	: P001, IBC02, R0
Transport category (RID)	: 2
Hazard identification number (RID)	: 33

#### 14.7. Maritime transport in bulk according to IMO instruments

Not applicable

## **SECTION 15: Regulatory information**

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

R001

#### 15.1.1. EU-Regulations

#### **REACH Annex XVII (Restriction List)**

Contains no substance(s) listed on REACH Annex XVII (Restriction Conditions)

#### **REACH Annex XIV (Authorisation List)**

Contains no substance(s) listed on REACH Annex XIV (Authorisation List)

#### **REACH Candidate List (SVHC)**

Contains no substance(s) listed on the REACH Candidate List

#### **PIC Regulation (Prior Informed Consent)**

Contains no substance(s) listed on the PIC list (Regulation EU 649/2012 concerning the export and import of hazardous chemicals)

#### POP Regulation (Persistent Organic Pollutants)

Contains no substance(s) listed on the POP list (Regulation EU 2019/1021 on persistent organic pollutants)

#### Ozone Regulation (1005/2009)

Contains no substance(s) listed on the Ozone Depletion list (Regulation EU 1005/2009 on substances that deplete the ozone layer)

#### VOC Directive (2004/42)

VOC content

: < 90 %

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

### Explosives Precursors Regulation (2019/1148)

Contains no substance(s) listed on the Explosives Precursors list (Regulation EU 2019/1148 on the marketing and use of explosives precursors)

#### Drug Precursors Regulation (273/2004)

Contains no substance(s) listed on the Drug Precursors list (Regulation EC 273/2004 on the manufacture and the placing on market of certain substances used in the illicit manufacture of narcotic drugs and psychotropic substances)

#### Seveso Directive (Disaster Risk Reduction)

Seveso III Part I (Categories of dangerous substances)	Qualifying quantity (tonnes)	
	Lower-tier	Upper-tier
P5c FLAMMABLE LIQUIDS Flammable liquids, Categories 2 or 3 not covered by P5a and P5b	5000	50000

#### 15.1.2. National regulations

No additional information available

#### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out

#### **SECTION 16: Other information**

Abbreviations and acronyms:		
ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways	
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road	
ATE	Acute Toxicity Estimate	
BCF	Bioconcentration factor	
BLV	Biological limit value	
BOD	Biochemical oxygen demand (BOD)	
COD	Chemical oxygen demand (COD)	
DMEL	Derived Minimal Effect level	
DNEL	Derived-No Effect Level	
EC-No.	European Community number	
EC50	Median effective concentration	
EN	European Standard	
IARC	International Agency for Research on Cancer	
ΙΑΤΑ	International Air Transport Association	
IMDG	International Maritime Dangerous Goods	
LC50	Median lethal concentration	
LD50	Median lethal dose	
LOAEL	Lowest Observed Adverse Effect Level	
NOAEC	No-Observed Adverse Effect Concentration	
NOAEL	No-Observed Adverse Effect Level	
NOEC	No-Observed Effect Concentration	
OECD	Organisation for Economic Co-operation and Development	
OEL	Occupational Exposure Limit	

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 SDS No: 00056-0072

Abbreviations and acronyms:		
РВТ	Persistent Bioaccumulative Toxic	
PNEC	Predicted No-Effect Concentration	
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail	
SDS	Safety Data Sheet	
STP	Sewage treatment plant	
ThOD	Theoretical oxygen demand (ThOD)	
TLM	Median Tolerance Limit	
VOC	Volatile Organic Compounds	
CAS-No.	Chemical Abstract Service number	
N.O.S.	Not Otherwise Specified	
vPvB	Very Persistent and Very Bioaccumulative	
ED	Endocrine disrupting properties	
DOT	Department of Transport	
TDG	Transportation of Dangerous Goods	
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006	
GHS	Globally Harmonized System of Classification, Labelling and Packaging of Chemicals	
IBC-Code	International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk	
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008	
MARPOL 73/78	MARPOL 73/78: International Convention for the Prevention of Pollution From Ships	
ADG	Transport of Australian Dangerous Goods	

Other information

\_\_\_\_

: Data of sections 4 to 8, as well as 10 to 12, do partly not refer to the use and the regular employing of the product (in this sense consult information on use and on product), but to liberation of major amounts in case of accidents and irregularities. The information describes exclusively the safety requirements for the product(s) and is based on the present level of our knowledge. The delivery specifications are contained in the corresponding product sheet. This data does not constitute a guarantee for the characteristics of the product(s) as defined by the legal warranty regulations.

Full text of H- and EUH-statements:		
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2	
Flam. Liq. 2	Flammable liquids, Category 2	
H225	Highly flammable liquid and vapour.	
H319	Causes serious eye irritation.	
H336	May cause drowsiness or dizziness.	
STOT SE 3	Specific target organ toxicity – Single exposure, Category 3, Narcosis	

Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:			
Flam. Liq. 2	H225	On basis of test data	
Eye Irrit. 2	H319	Calculation method	

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should therefore not be construed as guaranteeing any specific property of the product.