Vasco[®] OP Underglove Green. Safe. Sensitive.

Your Double Gloving System for a quick detection of perforations



CoE Practitioner Use





Vasco[®] OP Underglove Green. Safe. Sensitive.

Green colour for more visibility and safety

Easy double donning in combination with Vasco[®] OP Sensitive due to adjusted surface and fit

Reduced wall thickness for high tactile sensitivity

Double Gloving to increase personal protection

Wearing two pairs of gloves during surgery reduces the risk of blood-borne infections for healthcare professionals and their patients. Double gloving helps to quickly detect outer glove perforations and decreases the risk of inner glove perforations significantly. In case of liquid entrance through the outer glove, a dark spot becomes immediately visible when combining the green Vasco[®] OP Underglove and a second white Vasco[®] surgical glove.

The green colour allows a quick detection of perforations and gloves can be changed instantly.

size	sales unit	REF
6	40 pairs	6081200
6.5	40 pairs	6081218
7	40 pairs	6081226
7.5	40 pairs	6081234
8	40 pairs	6081242
8.5	40 pairs	6081259
STERILE	R (E 0123	

- Surgical glove, green
- Natural rubber latex
- Powder-free with synthetic inner coating
- Wall thickness at palm: min. 0.17 mm
- AQL 1.0



All Vasco® OP gloves are sterile, packed as right/left pairs in a paper wallet and peel pouch.